# InSpec<sup>™</sup> HA

Comprehensive Technical Resource Sheet for the Full Range of InSpec HA Liquid Formats



#### Features

- Non-hazardous yet fully effective sporicide
- Very low residue, corrosion and odour
- 6 months validated 'in-use' shelf life
- All formats are presented as Ready-to-Use (RTU)
- Multiple bags for cleanroom transfer
- Ideal rotational partner for InSpec QT, InSpec QT+ and InSpec AN
- Manufactured in accordance with GMP
- Active substance compliant with Article 95 of the Biocidal Products Regulation (BPR)

#### Formulation

InSpec HA is a sterile sporicidal preparation of active chlorine released from hypochlorous acid.

Active substance concentration, Active Chlorine 500 ppm.

#### **Instructions for Use**

InSpec HA is designed for spraying, wiping and mopping applications.

Spray Bottles: Hold approximately 15cm to 20cm from area to be treated. Apply to surface to ensure complete coverage for the required contact times.

Screw Cap Versions: Pour into an appropriate container for mopping. Apply to surface to ensure complete coverage for the required contact times.

#### **Material Compatibility**

Application of solutions, when used as directed, will not affect materials



normally encountered in the cleanroom. See compatibility information in the technical file.

#### Microbiological Minimum Efficacy Contact Times

Standard	Contact Time
EN 1276 Bacteria	5 Minutes
EN 13697 Bacteria	5 Minutes
EN 1650 Fungi	15 Minutes
EN 13697 Fungi	15 Minutes
EN 13704 Spores	1 Minute

#### Safe Handling and Storage Information

Always wear gloves and goggles or face protection. Always read the label and SDS before use.

Store upright in original closed containers, away from sunlight and extremes of temperature. Full guidance on the handling and disposal of this product is available in the Safety Data Sheet (SDS).

#### Formulation Batch Release Specifications

Specification	Release Parameters
SG	0.990 - 1.018
Concentration	350 - 500 ppm
рН	4.0 - 6.5
Colour	Colourless
Clarity	Clear
Odour	Faint chlorine odour

#### **Manufacturing Process**

InSpec HA is manufactured in accordance with GMP in an ISO 5 (Grade A) cleanroom. The solution is filtered through a 0.2-micron filter at point of fill.

#### Sterility

The sterile formats are aseptically filled into pre-gamma irradiated packaging (validated dose range of 25 - 45 kGy) to give a sterility assurance level (SAL) of 10<sup>-6</sup>.

Sterile formats of InSpec HA are tested for sterility for batch release.

The trigger spray formats have a validated 6-months 'in-use' sterility. This is delivered through a membrane-filter trigger.

For screw cap formats, use the entire contents in one session/4-hours to ensure sterility.

#### Certificates

InSpec HA is provided with a Certificate of Analysis (COA) confirming batch release specifications and providing batch manufacturing information.

The COAs for Sterile InSpec HA products contain information regarding the irradiation of the batch and the sterility test results.

### Technical Datasheet

## **In**Spec<sup>™</sup>**HA**



#### Formats Available

The following formats of InSpec HA are available:

Product	Code	Сар	Case Size	Container Material	Bags and Material	
Sterile Formats						
InSpec HA 1L Trigger Spray	HAPW50-1LS	Trigger	6 x 1L Bottles	HDPE	Double-Bagged LDPE	
InSpec HA 1L Dose	HAWFIDOSE-1LS	Screw Cap	6 x 1L Bottles	HDPE	Double-Bagged LDPE	
InSpec HA 5L Screw Cap	HAPW50-5LS	Screw Cap	2 x 5L Bottles	HDPE	Double-Bagged LDPE	
Non-Sterile Formats						
InSpec HA 1L Non-Sterile	NSHAPW-1L	Trigger	6 x 1L Bottles	HDPE	Double-Bagged LDPE	
InSpec HA 5L Non-Sterile	NSHAPW-5L	Screw Cap	2 x 5L Bottles	HDPE	Double-Bagged LDPE	

#### **Transport Information**

Product	Code	Case Dimensions	Commodity Code	Cases per Pallet Euro/UK	Dangerous Goods	
Sterile Formats	Sterile Formats					
InSpec HA 1L Trigger Spray	HAPW50-1LS	27cm x 19cm x 33cm 6.4Kg	29159070	64/60	Not classified as hazardous for transport	
InSpec HA 1L Dose	HAWFIDOSE-1LS	27cm x 19cm x 33cm 6.4Kg	29159070	64/60	Not classified as hazardous for transport	
InSpec HA 5L Screw Cap	HAPW50-5LS	29cm x 20.5cm x 31cm 10.8Kg	29159070	44/51	Not classified as hazardous for transport	
Non-Sterile Formats						
InSpec HA 1L Non-Sterile	NSHAPW-1L	27cm x 19cm x 33cm 6.4Kg	29159070	64/60	Not classified as hazardous for transport	
InSpec HA 5L Non-Sterile	NSHAPW-5L	29cm x 20.5cm x 31cm 10.8Kg	29159070	44/51	Not classified as hazardous for transport	



Discovery 2 • NETPark • 2 William Armstrong Way • Sedgefield County Durham • TS21 3FD • UK UK Tel: 0845 073 9529 • International Tel: +44(0) 1527 830 942 Web: www.redditchmedical.com • Email: sales@redditchmedical.com HA Group Tec